

AH Plus® Bioceramic Sealer

**Ready, SET, Go!**




Dentsply  
Sirona

#EndoPassion

AH Plus® Bioceramic Sealer

## What sets this sealer apart?

Bioactivity of bioceramic sealers is important but you need more to ensure a stable seal. The apical seal is only completed after the sealer sets. AH Plus Bioceramic sealer sets apart from the rest because it sets quickly, has a high washout resistance and is more radiopaque compared to current market leaders.



1. data on file, compared to EndoSequence BC Sealer, Total Fill BC Sealer, EdgeBioceramic Sealer



- 1 A faster set**  

4 hours Our Bioceramic Sealer  
10 hours EndoSequence BC Sealer™  
Total Fill BC Sealer™
- 2 Dimensionally stable**  

Low solubility of 0.11% with a high washout resistance to stably seal the root canal for a long time.

Solubility  
AH Plus BC 0.11%  
EndoSequence BC Sealer 3% up to 23%  
Total Fill BC Sealer

1. data on file. Same solubility results for EndoSequence BC Sealer, Total Fill BC Sealer, EdgeBioceramic Sealer.
- 3 Safe and biocompatible**  

Induces tissue healing and does not stimulate periodontal tissue, even when in contact.<sup>2</sup>
- 4 Better radiopacity**  

25% more radiopaque than the EndoSequence and TotalFill BC Sealer.<sup>2</sup>

AH Plus Bioceramic Sealer EndoSequence BC Sealer

2. data on file.
- 5 Enhanced aesthetics**  

Free of Bismuth Oxide, our sealer does not discolor the tooth and guarantees a confident smile.<sup>3</sup>
- 6 Retreatability**  

Removable even after setting with a general hand file or NiTi file.

3. formulation of product

## AH Plus® Bioceramic Sealer



Starter Kit	B KT AHBS 000 000	Syringe Refill	B 00 AHBS 100 03G	Cannula Refill	B 00 AHBS 020 24G
-------------	-------------------	----------------	-------------------	----------------	-------------------

### ProTaper Ultimate™ Conform Fit® Gutta-Percha Points

#### Box of 60 gutta-percha points

<b>F1</b>	B NS PULG FGP OF1
<b>F2</b>	B NS PULG FGP OF2
<b>F3</b>	B NS PULG FGP OF3
<b>FX</b>	B NS PULG FGP OFX
<b>FXL</b>	B NS PULG FGP FXL
<b>ASSORTMENT</b>	B NS PULG FGP AST
<b>F1-F3</b>	B NS PULG FGP AST
12x F1, 24x F2, 24x F3	

### WaveOne® Gold Conform Fit™ Gutta-Percha

#### Box of 60 variable gutta-percha

<b>Small</b>	B 00 WGGP F00 SML
<b>Primary</b>	B 00 WGGP F00 PRI
<b>Medium</b>	B 00 WGGP F00 MED
<b>Large</b>	B 00 WGGP F00 LRG
<b>ASSORTMENT</b>	B 00 WGGP F00 AST


### TruNatomy™ Conform Fit™ Gutta-Percha Obturation Points

#### Box of 60 gutta-percha points

<b>Small</b>	B 00 TNGP F00 SML
<b>Prime</b>	B 00 TNGP F00 PRM
<b>Medium</b>	B 00 TNGP F00 MED
<b>ASSORTMENT</b>	B 00 TNGP F00 AST

## Triple lock your obturation

Only DS treatment solutions can give you a unique Triple Lock™: Ensure an ideal shape for your workflow. Get an unrivaled fit with Conform Fit® Gutta-Percha points matching these unique shapes.¹ Finally, achieve a gap-free bond to gutta-percha and the canal walls with AH Plus® Bioceramic Sealer.²



### Bioceramic Seal

AH Plus® Bioceramic Sealer quickly sets, the final step towards inducing tissue healing.³

### Precise Shape

Every tooth is different; every clinician has preferences. The files of our treatment solutions ensure an ideal shape for your workflow.

### Conform Fit® Technology

Less waste and saves time because the first cone fits every time.¹

## Discover which DS treatment solution fits best.

ProTaper Ultimate™ | WaveOne® Gold | TruNatomy®

<b>Gutta-Percha</b>	<b>Dentsply Tulsa Dental Specialties</b>	<b>DE</b>	<b>CE</b>
	808 Rolling Hills Drive, Johnson City, TN 37604	<b>Dentsply DeTrey GmbH</b>	2757
<b>AH Plus® Bioceramic Sealer</b>	Harsanki 2/208, Medical Industry Complex Strada - 43-10, Tanjungpinang-gil, Wörnu-si, Ganganon-da, KR-26311	<b>Dentsply DeTrey GmbH</b>	0937
		<b>Biogel</b>	Am Pflanzhof 12, DE-22303 Hamburg



INL // BEN AHPO BRC 000 / 00/2/2021 - created 07/2021

## Set your new standard

### Ah Plus® Bioceramic Sealer

- 1 Biocompatible**  
Providing an ideal environment for Hydroxyapatite formation paving the way for the body's self healing process.¹
- 2 Gap-free**  
Our sealer produces a wide range of hydroxyapatite crystals within the dentinal tubules. Being lipophilic and hydrophilic, leading to a gap-free bond with the canal wall and gutta-percha simultaneously. No need for surface treated gutta-percha points.¹
- 3 Wash-out resistant**  
Thanks to its low solubility and low film thickness our sealer stably seals the root canal.
 

Film Thickness²		Solubility³	
AH Plus BC	25	AH Plus BC	0.1%
EndoSequence BC Sealer	75	EndoSequence BC Sealer	3% upto 23%
- 4 Convenience**  
Convenience delivered in a pre-mixed syringe with a small 24 gauge that can be placed directly into the canal.

### Conform Fit® Technology

100% fit⁴

Range of fit

vs

Range of fit with Conform Fit® Gutta-Percha Point

21% fit⁵  
When using Cheap Endo gutta-percha cones

Range of fit

vs

Range of fit with Cheap Endo gutta-percha

1. data on file 4. on 28 mold cavities 5. 28 random cones